





كیر میدیکال تریدینگ ذ.م.م  
CARE MEDICAL TRADING W.L.L



# IN CONTINENCE AND UROLOGY PRODUCTS



كیر میدیکال تریدینگ ذ.م.م.  
CARE MEDICAL TRADING W.L.L



# FOI FY CATHETERS



## Description

Size : 6,8,10Fr/Ch  
Balloon Capacity : 6Fr/Ch 3cc,8,10Fr/Ch 3-5cc  
Length : 250mm  
Made from natural latex. Silicone Coated.  
Available with the option of plastic valve or rubber valve.  
Natural or Hydrophilic.  
Packed sterile.For Single use only.



## Description

Size : 12,14,16,18,20,22,24,26,28,30Fr/Ch  
Balloon Capacity : 12-24Fr/Ch 5-10cc,5-15cc  
: 12-30Fr/Ch 30cc  
: 16-30Fr/Ch 30-50cc  
Length: 400mm(Male) or 270mm(Female)  
Made from natural latex. Silicone Coated.  
Natural or Hydrophilic and Male or Female.  
Packed sterile. For Single use only.



## Description


Size : 16,18,20,22,24,26Fr/Ch  
Balloon Capacity : 30cc  
Length : 400mm  
Made from natural latex. Silicone Coated.  
Natural or Hydrophilic.  
Packed sterile. For Single use only.

# SILICONE FOLEY CATHETERS



## Description


Size : 8,10,12,14,16,18,20,22,24Fr/Ch  
Balloon Capacity : 8-10Fr/Ch 3cc,3-5cc,5cc  
: 12-24Fr/Ch 5-10cc,5-15cc,5cc,10cc,15cc  
: 16-24Fr/Ch 30cc  
Length : 400mm  
Made from 100% Silicone.  
Packed sterile. For Single use only



## Description

Size : 16,18,20,22,24Fr/Ch  
Balloon Capacity : 5-15cc,30cc  
Length : 400mm  
Made from 100% Silicone.  
Packed sterile. For Single use only.

# SUCTION CATHETER





كیر میدیکال ٹریڈینگ ذ.م.م  
CARE MEDICAL TRADING W.L.L

TM  
greaturo

## NELATON CATHETER



### Description

Size : 6,8,10,12,14,16,18,20,22,24Fr/Ch  
Length : 400mm(Male) or 180mm(Female)  
Made from PVC or SOFT PVC.  
Natural or Hydrophilic and Male or Female.  
Packed sterile. For Single use only.

## ECONOMY URINE BAG

### Description



Ref : GCU201501  
Front Film : Transparence  
Back Film : White  
Outlet : Pull-push valve  
Capacity : 2000ml  
A.1500ml and 1000ml are available.  
B. With anti-reflux valve. Tubing lengh will be 90cm or 110cm.  
C. Available in two plastic hangers or belts. Different distance is available between double holes.  
D. Latex-free.Packed in poly bag or blister with E.O. sterile. Or bulk. Greatcare brand or OEM available.



Ref : GCU201505 With reinforced plastic hang tape  
Front Film : Transparence  
Back Film : White  
Outlet : Pull-push Valve  
Capacity : 2000ml



Ref : GCU201701  
Front Film : Transparence  
Back Film : White  
Outlet : Without  
Anti-reflux valve: Without  
Capacity : 2000ml

Ref : GCU201703  
Front Film : Transparence  
Back Film : White  
Outlet : Without  
Anti-reflux valve: With  
Capacity : 2000ml

A.1500ml and 1000ml are available.  
B. With anti-reflux valve. Tubing lengh will be 90cm or 110cm.  
C. Available in two plastic hangers or belts. Different distance is available between double holes.  
D. Latex-free.Packed in poly bag or blister with E.O. sterile.  
Or ulk. Greatcare brand or OEM available.



كير ميديكال تریدينغ ذ.م.م  
CARE MEDICAL TRADING W.L.L

greauro™

## DELUXE URINE BAG

### Description

Urine Bag  
2000ml,with white sample port & T valve.  
Extra thickness film(0.15mm per ply),double transparent films.  
Extra outer diameter tubing(9mm).  
Four holes for different distance hang sizes.



### Description

Urinary Drain Bag  
2000ml,with T valve.  
Front film is transparent, back film is white.  
With white sample port.With anti-reflux valve.  
Latex-free.Packed sterile.



### Description

Urinary Drain Bag  
2000ml,with twist valve.  
Front film is transparent, back film is white.  
With sample port. With anti-reflux valve. With plastic tie.  
Latex-free.Packed sterile.



### Description

Urinary Drain Bag  
2000ml,with T valve.  
Front film is transparent, back film is white.  
With sample port.With anti-reflux valve.With plastic tie.  
Latex-free.Packed sterile.



### Description

Urinary Drain Bag  
2000ml,with twist valve.  
Front film is transparent, back film is white.  
With sample port.With anti-reflux valve. With reinforced single hanger.  
Latex-free.Packed sterile.





كیر میدیکال تریدینگ ذ.م.م.  
CARE MEDICAL TRADING W.L.L

greaturo™

# DELUXE URINE BAG



## Description

Urinary Drain Bag  
2000ml,with twist valve.  
Front and back film are transparent.  
With anti-reflux valve.With plasic tie.  
Latex-free.Packed sterile.



## Description

Urinary Drain Bag  
2000ml,with T valve. With white sample port.  
Front film is transparent, and back film is white.  
With anti-reflux valve. With reinforced double hanger. With  
tubing clamp.Elastic tubing.  
Latex-free. Packed sterile



## Description

Urinary Drain Bag  
Ref.: GCU210701 2600ml  
Ref.: GCU210702 2500ml  
Ref.: GCU210703 2000ml  
With anti-reflux tower(three part),With needleshield sample port.  
With bed sheet clamp. Tubing clip is optional.  
With reinforced double hanger & rope hanger.  
With T valve.Latex-free.Packed sterile.  
Blister pack or poly bag packing.



## Description

Urinary Drain Bag  
Ref.: GCU210801 2600ml  
Ref.: GCU210802 2500ml  
Ref.: GCU210803 2000ml  
With anti-reflux tower(one part),With needleshield sample port.  
With bed sheet clamp. Tubing clip is optional.  
With reinforced double hanger & rope hanger.  
With T valve.Latex-free.Packed sterile.  
Blister pack or poly bag packing.



## Description

Urinary Drain Bag  
2000ml, With anti-reflux tower, With needleshield sample port.  
With reinforced single hanger. With finger guard clamp.  
Latex-free.Packed sterile.  
Blister pack or poly bag packing.



كیر میدیکال تریدینگ ذ.م.م  
CARE MEDICAL TRADING W.L.L

greauro™

## DELUXE URINE BAG

### Description



Urinary Drain Bag  
2000ml, With anti-reflux tower, With white sample port.  
With reinforced single hanger. With T valve outlet.  
Latex-free.Packed sterile.  
Blister pack or poly bag packing.

### Description



Urinary Drain Bag  
Ref: GCU211701 2600ml  
Ref: GCU211702 2500ml  
Ref: GCU211703 2000ml  
With anti-reflux device(two parts), with needless sample port.  
With bed sheet clamp. Tubing clip is optional.  
With reinforced double hanger & rope hanger.  
With T valve.Latex-free.Packed sterile.  
Blister pack or poly bag packing.

## URINE METER



Ref.: GCU220101 2600ml drain bag and 200ml urine meter.  
Ref.: GCU220102 2500ml drain bag and 200ml urine meter.  
Ref.: GCU220103 2600ml drain bag and 400ml urine meter.  
A. Easy-to-read molded plastic meters are calibrated up to 200ml.  
B. With anti-reflux valve.With needless sample port.  
C. With reinforced double hanger and rope hanger.  
D. Latex-free.Packed sterile.



2600ml drain bag and 500ml urine meter.



2600ml drain bag and 400ml urine meter.  
With big reinforced double hanger & rope hange.



كیر میدیکال تریدینگ ذ.م.م  
CARE MEDICAL TRADING W.L.L

## DELUXE URINE BAG




2600ml drain bag and 500ml urine meter.  
With drip chamber and white sample port.



2000ml drain bag and 200ml urine meter.

greaturo™

## LEG BAG



500ml, ellipse shape, with  
green twist valve &  
green tip cover.  
With double high elastic belts.



600ml, square shape, with  
green pull-push valve &  
green tip cover. With double high  
elastic belts.



500ml, one ply flocking fabric,  
another ply is PVC film.  
With T valve & short tubing.



500ml, one ply flocking fabric,  
another ply is PVC film.  
With T valve & long tubing.




500ml, one ply flocking fabric,  
another ply is PVC film.  
With double reversal valve &  
long tubing.




كير ميديكال تریدينغ ذ.م.م  
CARE MEDICAL TRADING W.L.L

greaturo™

## URINE COLLECTOR



Pediatric Urine Collector  
100ml, frankseal,  
yellow paper.



Pediatric Urine Collector  
100ml, frontseal, yellow  
paper.

## ENEMA SET



Pediatric Urine Collector  
160ml, frankseal, white  
paper.



Enema Set 1500ml

## WOUND DRAINAGE RESERVOIR



Wound Drainage Reservoir  
With anti-reflux valve, with different sizes needle,  
such as 8fr, 10fr, 12fr, 14fr, 16fr, 18fr, etc.





كیر میدیکال تریدینگ ذ.م.م  
CARE MEDICAL TRADING W.L.L

# URETHRAL CATHETER



Size: ,8,10,12,14,16,18,20,22,24,26,28Fr/Ch  
Length: 400mm  
Packed sterile.For Single use only.



Size: 6,8,10,12,14,16,18,20,22,24,26,28Fr/Ch  
Length: 400mm  
Packed sterile.For Single use only.

greaturo™

# MALECOT CATHETER



Size: 12,14,16,18,20,22,24,26,28,30,  
32,34,36,38,40Fr/Ch  
Length: 400mm  
Packed sterile.For Single use only.

# MALE EXTERNAL CATHETER



Packed sterile, For single use only



100% silicone,self-adhering,  
with white tray.  
transparent tray available



100% Silicone,Self-adhering  
Type: Standard,Pop on, Wide Band  
Size: 24mm,28mm,31mm,35mm,40mm



كیر میدیکال تریدینگ ذ.م.م  
CARE MEDICAL TRADING W.L.L

## URINE BAG HANGER



## URINE BAG VALVE



greauro™

## URINE BAG CLAMP



## URINE BAG CONNECTOR



## DIAPER



DIAPER X-LARGE WITH INDICATOR  
DIAPER LARGE WITH INDICATOR, 80/CTN  
DIAPER SMALL WITH INDICATOR, 80/CTN  
DIAPER MEDIUM WITH INDICATOR, 80/CTN  
DIAPER (PAMPERS) 16+  
DIAPER (PAMPERS) 9-20KG  
DIAPER (PAMPERS) 4-9 KG  
DIAPER (PAMPERS) 2-5 KG  
DIAPER NEW BORN 3-6 KG 64 X 3/CTN  
DIAPER(PAMPERS) 6X18/CTN



## UNDERPAD

UNDER PAD 60 X 45 CM 38G  
UNDER PAD 60 X 60 CM 60G-240/CTN  
UNDER PAD 60 X 90 CM70G-150/CTN  
UNDER PAD 60 X 90 CM80G-150/CTN



كير ميديكال تریدينغ ذ.م.م  
CARE MEDICAL TRADING W.L.L

## CATHETERIZATION TRAY

greaturo<sup>TM</sup>



كير ميديكال تریدينغ ذ.م.م  
CARE MEDICAL TRADING W.L.L

Care Medical Trading, Office No:205, Bin Dirham Plaza  
Street 220, B-Ring Road, Zone 25, Al mansoura, Doha-Qatar  
+974 44120444 / sales@caremedicals.com / www.caremedicals.com