



كير ميديكال تریدينغ ذ.م.م
CARE MEDICAL TRADING W.L.L



MEDICAL FURNITURES

EXAMINATION COUCH AND BED



Stainless steel frame
Mechanical back adjustment (optional)
Feet with rubber cushion
Upholstered with polyurethane rigid foam covered with washable PU material, easy to clean, hard wearing

EXAMINATION COUCH AND BED



Examination table
1.CE&ISO13485 Examination table
2.Powder coated steel examination table
3.Hydraulic foot pedal operation



EMERGENCY CART

Description: Emergency Trolley

High-strength ABS boards ,high-quality metal ball-bearing rail.

Telescopic IV pole,

Defibrillator shelf,

Heart board,

Power outlet&Hooks,

Oxygen tank holder

Waste Containers ,

Utility container,

File bag, sliding

side shelf,

sharps container holder.

Top with stainless steel guard rail Aluminum columns with impact protection design

Centralized Locking system

Work surface top covered with transparent soft plastic glass

3 big & 2 small, inner with dividers can be organized freely into different size units

4 anti-static wheels (Dia:125MM), two with brakes



MEDICAL CABINET



BED SIDE CABINET



IV STAND

Available in full stainless steel/ epoxy powder coated stand & stainless steel saline rod.

